YON DON[®]

NFRARED WARMING LIGHT

The **YD-WL-DTL** is aimed to provide professional therapy needs which were not available until Based on YON DON's expertise and 26 years of experience in manufacturing medical now. equipment, the **YD-WL-DTL** features safety, durability, ease of use, and especially the extraordinary lamp arm which allows therapy application from high, low, far, and close-reach from the lamp body.



Figure 1: Common application

Specifications:

- 1. Voltage: 100V ~ 110V @ 50/60 Hz. 220V ~ 240V @ 50/60 Hz.
- 2. Radiating power: 0 ~ 175W, 0 ~ 250W, 0 ~ 375W.
- 3. Mobile stand with castors fixed at base.
- 4. Reinforced lamp socket for durability and safety.
- 5. Good heat dispatching lampshade with stainless steel protective screen, which prevents burning accidents.
- 6. Adjustable arm height: 60 cm to 150 cm, and adjustable lamp post height.
- 7. Adjustable lamp arm for various extension styles (Figure 3).
- 8. At 80 cm from the lamp post, the therapy unit still remains stable and firm.
- 9. Adjustable radiating angle within 360 degrees.
- 10. Both timed operation and continuous operation are available for operators' convenience. Built-in auto-sleep timer allows timed operations for up to 60 minutes (or longer, upon customer's request).
- 11. Built-in dimmer for adjusting radiating power.
- 12. Can be easily retracted to a more compact size for space saving storage and/or close-reach application (Figure 2).
- * The above specification is subject to changes without notice.

Manufacturer, Exporter & Importer YON DON ENTERPRISE CO., LTD.

2F, NO. 34, LANE 80, NAN KANG RD., SEC. 3, TAIPEI 11570, TAIWAN. R.O.C. TEL: 886-2-2783-9583, 886-2-2785-9800 (Rep.) FAX: 886-2-2786-5539

http: www.yondon.com e-mail: info@vondon.com



Figure 2: Retracted for space-saving storage and for close-reach application.



Figure 3: Therapy heating source may be moved to lower position for applications at or below knee point.